

# EasyOne Air

*Portable and PC spirometer*



EasyOne Air is a flexible, accurate and proven spirometry solution. Easily connect to the PC with Bluetooth for real-time patient incentives, data exchange and EMR connectivity.

# Simple. Certain. Proven.



## Features

- Large color touch screen
- Easy navigation
- Rechargeable battery
- Portable or PC mode
- Real-time animated incentives
- Bluetooth connectivity
- Flexible use, consistent results
- Real-time curve display
- Real-time data transfer
- Proven EMR integration with top EMR/EHR systems
- Calibration free technology



## Specifications

**Test Types:** FVC, FVL, SVC, MVV and Pre/Post Testing

**Languages:** English, French, German, Spanish, Swedish

**EMR connectivity:** EasyOne Connect

**Dimensions:** 3.4" x 6.1" x 1.4", 13 oz

**Standards:** Meets ATS/ERS, NIOSH/OSHA and SSA Disability

## Part Numbers

**2500-2A:** EasyOne Air spirometry system

**5050-50:** EasyOne Flow Tube box of 50

**5050-200:** EasyOne Flow Tube box of 200

  
**TrueFlow**

The original ultrasonic flow measurement is highly accurate in all flow ranges, independent of gas composition, pressure, temperature and humidity.

n d d TrueFlow is a hygienic and resistance-free solution that does not require calibration during its lifetime.

  
**EasyOneConnect**

EasyOne Connect software offers one database and one platform for all EasyOne point-of-care-solutions. Data management has never been easier.

