



## Amplivox Otowave Tympanometer Comparison Matrix

User Features	102/102-4	202/202-H	302/302-H
Portable	✓	✓	✓
PC connectivity (USB)	-	✓	✓
PC connectivity (IrDA)	✓	-	-
Internal memory	✓	✓	✓
NOAH compatible	✓	✓	✓
Thermal printer output	✓	✓	✓
Programmable 1 frequency (1kHz) ipsilateral reflex	✓ / -	-	-
Programmable 4 frequency ipsilateral reflex function	- / ✓	✓	✓
Programmable 4 frequency contralateral reflex function	-	✓	✓
Ergonomic design	✓	✓	✓
Mains power function (adaptor supplied)	-	✓	✓
Battery power function	✓	✓	-
Carry case supplied	✓	✓	✓
Daily test cavity pack supplied	✓	✓	✓



## Amplivox Otowave Tympanometer Comparison Matrix

Technical Features	102/102-4	202/202-H	302/302-H
226Hz probe tone +/-0.1%	✓	✓	✓
1000Hz probe tone +/-0.1%	-	- / ✓	- / ✓
User selectable sweep speeds	-	✓	✓
+200daPa to -400daPa pressure range +/- 10daPa	✓	✓	✓
0.2 to 5ml volumetric range +/-0.01ml	✓	✓ / -	✓ / -
0.1 to 5ml volumetric range +/-0.01ml	-	- / ✓	- / ✓
1kHz ipsi reflex function	✓ / -	-	-
500, 1k, 2k and 4k ipsi reflex function	- / ✓	✓	✓
500, 1k, 2k and 4k contralateral reflex function	-	✓	✓
70 to 100dB reflex range	✓	✓	✓
0.01mLl to 0.5ml reflex measurement range	✓	✓	✓
Weight (grams)	380 g	430 g	760 g
Dimensions wide x depth x height (mm)	80 x 210 x 40	85 x 190 x 40	270 x 175 x 70