

ER-75

OUTPUT DATA

| | |
|-----------------------|-----------------|
| Frequency range (Hz): | 250Hz to 8000Hz |
| Threshold (1KHz): | 75dBHL |
| Repetition: | 3dB |
| Hold time: | 330ms |

PHYSICAL DATA

| | |
|------------------------|------------------------------|
| Power (battery): | 2 x 1.5V "AA" type batteries |
| Construction: | ABS plastic enclosure |
| Operating temperature: | 0°C to 40°C |
| Dimensions (mm): | 158 x 110 x 115 high |
| Weight: | 500g |

EQUIPMENT

STANDARD

- Operating manual
- Interconnecting lead
- 2 x AA batteries



For more information, please contact

e: sales@kahntactmedical.com | t: 800.742.4909 | w: kahntactmedical.com

The Amplivox policy is one of continuous development and consequently the equipment may vary in detail from the description and specification in this publication. Product specification may differ by country.