


**Alcovisor
Mark V**



**ALCOHOL TESTING
INSTRUMENTS**

- One Year Factory Warranty
- Perform a test on a matter of seconds
- Low-cost
- Store 10,000 test records in memory
- External thermal printer
- NHTSA approved Evidential Breath Testing Device

800.742.4909
kahntactmedical.com

Alcovisor Mark V Breathalyzer

Quick DOT Testing

The Mark V is approved by NHTSA as an Evidential Breath Tester and approved for DOT workplace testing. Fast and easy to use, the Mark V conducts a breath alcohol test in less than 60 seconds.

Mark V Product Benefits & Features

- Able to store 10,000 test records in its memory
- Fast testing in a matter of seconds
- Affordable
- Automatic and manual breath sampling enabled
- Color LCD screen
- Evidential accuracy
- External thermal printer
- Audible beeping as subject provides a breath sample let's you know air is flowing through the device
- Visual bar on the LCD screen shows you how much breath the subject is providing
- Automatic Air Blank taken at the start of each test
- USB PC Connection available with optional upload software allows you to upload test records to a computer

Stores Thousands of Test Results

The robust memory of the Alcovisor Mark V allows the instrument to store up to 10,000 test results. The operator can browse the test results stored in the instrument and re-print selected tests, or upload results to a computer with the optional upload software.



RELATED ITEMS AND SUPPLIES

Item	Description	Price
389-MARK	Mouthpiece for the Mark V	\$.31 @ 1-499
		\$.28 @ 500-999
		\$.25 @ 1,000-4,999
388-THERM	Alcovisor Printer Paper (4 rolls)	\$13.50
324-KT	Alcohol Gas Tank, 105 liter	\$170.00
382-04KT	Alcohol Gas Tank, 34 liters	\$123.00
328	Regulator	\$112.00
371	DOT Alcohol Test Forms	\$0.25 ea.
372	Non-DOT Alcohol Test Forms	\$0.25 ea.
342-KT	Yellow Tamper Tape (1 box)	\$25.00

STANDARD PRODUCT

Item	Description	Price
200	Mark V with Thermal Printer Includes: Device, printer, standard carrying case, 105 mouthpieces, batteries, printer paper, passive sampling cup, manuals, 1-yr warranty	\$985.00
200-KIT	Mark V with Thermal Printer (DOT Package) Includes: Device, printer, large carrying case, 105 mouthpieces, batteries, printer paper, passive sampling cups, manuals, 34-liter gas tank, 100 DOT and non-DOT alcohol testing forms, box of tamper tape, 1-yr warranty	\$1,285.00

SPECIFICATIONS

Sensor type: Electrochemical Oxidation (fuel cell)

Range of measurement: .000 - .400 BrAC

Accuracy and precision: U.S. DOT approved for evidential use

Breath sampling: Automatic with manual override

Response time: 3-5 seconds negative results; 5-10 seconds positive results

Power requirement: AAA batteries (included) powers device. Lithium battery powers the printer (included)

Computer capability: Optional upload software

Calibration: Microprocessor controlled, using gas or wet bath simulator

Memory: 10,000 tests

Print capability: With external thermal printer

Working temperature range: -10 to 50°C

Warranty: One year parts and labor

For Information & Orders Call: **800.742.4909**
Order 24 Hours A Day At: www.kahntactmedical.com
Fax Orders To: **708.449.7832**

Quality Medical Devices Since 1987
Knoxville, TN 37939