

Vitalograph Spirometer Selection Guide	PC Spirometers		Handheld/Portable Spirometers			Desktop/Portable
	Pneumotrac Spirotrac 6 Software	Pneumotrac Spirotrac V Software	Micro Reports Software	In2itive Reports Software	In2itive Spirotrac V Software	ALPHA Spirotrac 6 Software
Meets 2019 ATS/ERS Spirometry Standards	Yes	No	Yes	Yes	No	Yes
Tests:						
VC	•	•	•	•	•	•
FVC	•	•	•	•	•	•
Post BD	•	•	•	•	•	•
MVV	•	•			•	•
Challenge Test	•	•			•	•
6MWT/COPD Assess	•	•			•	•
Total Parameters Available	50	60	8	50+	60	60
Quality Measures and Incentives:						
Predicted Values	•	•	•	•	•	•
Interpretation	•	•		•	•	•
Quality Prompts	•	•	•	•	•	•
Incentives	•	•		•	•	•
Data Exchange and Printing:						
Built-in Printer						•
PDF Reports (S = Spirotrac; V = Reports software; DS= Device Studio software)	•S	•S	•V	•V	•S/V	•DS
Bi-Directional Data Exchange with EMRs	•	•			•	•
General Features:						
Touch-Screen			•	•	•	•
Color Display		•	•	•	•	•
Internal Battery			•	•	•	•
Full Database	•	•		•	•	•
Add-On (Testing Options):						
Pulse Oximetry	•	•			•	
Resting 12-Lead ECG	•	•			•	